

## Culture flasks on slide



### Nunc™ culture flasks on slide



Flaskette (white cap)

Flask (black cap)

#### Applications

- Autoradiography
- Single-cell immunofluorescence

Cat. No.	Description	Packing units	Units/ carton	€ Excl. VAT/ carton
055065	Flask	5	50	NC -
055066	Flaskette	8	16	NC -

#### Flask

- ▶ Flask attached to a Nunclon™ Delta culture-treated polystyrene slide by ultrasonic welding for long-term, leak-free incubation
- ▶ The flask can be easily detached from the slide

#### ▶ Recommended volume depending on the surface

Culture area (cm <sup>2</sup> )	Working volume (ml)/chamber
10	2,5- 5

#### Flaskette

- ▶ Flask (20 x 50 mm) fixed on a glass slide with a seal
- ▶ The flask can be detached from the slide

#### ▶ Recommended volume depending on the surface

Culture area (cm <sup>2</sup> )	Working volume (ml)/chamber
9	2,5- 5



### Hybridwell™ culture flasks on SPL slides



#### ▶ Sterile

- ▶ Polystyrene chamber glued to a glass microscope slide (75 x 25 mm)

#### ▶ Non-cytotoxic

#### ▶ Slide available in different materials:

- glass, naturally charged, ideal growth surface for adherent
- DLux, modified plastic with minimal autofluorescence ;
- FLux, modified plastic film (188 µm) for confocal microscopy



#### ▶ Recommended volume depending on the surface

Culture area (cm <sup>2</sup> )	Working volume (ml) / chamber
9	2.5-5,5

Cat. No.	Culture-treated	Packing units	Units/ carton	€ Excl. VAT/ carton
330069	No	6	12	NC -
330070	Yes	6	12	NC -
330071	Yes	6	12	NC -