

## Omni stainless steel dispersion rods with rotor stator rod

### Rods for Omni TH and GLH 850 homogenisers

Cat. No.	Dimensions Ø x h (mm)	Volume (ml)	Bottom	Type of grinding	€ Excl. VAT
078007	5 x 75	0.02 - 5	Flat	Fine	NC -
078008	5 x 75	0.02 - 5	Flat	Coarse	NC -
078148	7 x 115	0.25 - 30	Serrated	Fine	NC -
078149	7 x 115	0.25 - 30	Serrated	Coarse	NC -
078150	10 x 115	1.5 - 100	Serrated	Fine	NC -
078151	10 x 115	1.5 - 100	Serrated	Coarse	NC -

### Rods for Omni GLH 850 homogeniser only

Cat. No.	Dimensions Ø x h (mm)	Volume (ml)	Bottom	Type of grinding	€ Excl. VAT
078021	20 x 195	100 - 2000	Serrated	Fine	NC -
078022	20 x 195	100 - 2000	Serrated	Coarse	NC -
078152	30 x 195	200 - 5000	Crenellated	Fine	NC -
078153	30 x 195	200 - 5000	Crenellated	Coarse	NC -

### Rods for THQ Omni homogeniser

Cat. No.	Dimensions Ø x H (mm)	Volume (ml)	Bottom	Units/ box	€ Excl. VAT/ box
078212	5 x 75	0.2 - 5	Serrated	1	NC -
078213	5 x 75	0.2 - 5		6	NC -
078214	7 x 100	0.25 - 30		1	NC -
078215	7 x 100	0.25 - 30		6	NC -
078216	10 x 110	1.5 - 100		1	NC -
078217	10 x 110	1.5 - 100		6	NC -

### Accessoire

Cat. No.	Description	€ Excl. VAT
078026	Dismantling tools	NC -



The photos on the right are not to scale with each other.

## Plastic dispersion rods : Omni Tip™



- Reusable and autoclavable up to 7 times
- To eliminate the risk of cross-contamination between samples
- Ideal for PCR, DNA or RNA blotting analysis
- The clear plastic rod tube shows that there is no sample left
- Can be cleaned, disinfected and sterilised for repeated use
- Two-piece construction for easy disassembly and cleaning
- 80% less protein residue compared to a stainless steel tool
- Compatible with all Omni homogenisers with EZ connection (GLH, TGH, THQ)
- Available in sterile version: avoids cross-contamination and transfer between samples
- Diameter 7 x 110 mm
- Recommended workload : 0.25 ml to 30 ml

### Two models:

- rods for standard dispersion of liquid mixtures, soft tissues or tissues requiring gentle homogenisation;
- rods for the disintegration of hard or frozen tissues.

### Adapter

Cat. No.	Description	Units/ box	€ Excl. VAT/ box
078042	Adapter for homogenisers Omni EZC and Fisher PowerGen	1	NC -



### Dispersion rods

In boxes	25	50	100	500
<b>Omni Tip™ standard 110 mm dispersion rods</b>				
Cat. No.	078031	078032	078033	078034
€ Excl. VAT	NC -	NC -	NC -	NC -
<b>Omni Tip™ Dispersion Rods for hard or frozen tissues 110 mm</b>				
Cat. No.	078035	078036	078037	078038
€ Excl. VAT	NC -	NC -	NC -	NC -
<b>Omni Tip™ Dispersion Rods for hard or frozen tissues 110 mm sterile</b>				
Cat. No.	078039	-	-	-
€ Excl. VAT	NC -	-	-	-