Scope single tube scanner



- ▶1D and 2D code scanner
- No software required, just plug and play
- Visualisation of the code in the file of your choice (Word, Excel, OneNote...)
- > Works with PC, tablet, iPad, smartphone
- > 2 versions: USB cable or Bluetooth
- >Available as a Bluetooth version with an Lenovo 10"Android tablet
- Dimensions W x D x H : 59 x 150 x 38 mm
- >Weight : 254 g

Cat. No.	Description	€ Excl. VAT
783959	USB Scope Reader	NC -
783980	Bluetooth Scope Reader	NC -
783981	Bluetooth Scope Reader with Lenovo Android Tablet	NC -

Perception Scanner for tubes on rack

- Ultra fast: read a 96 rack in less than a second
- Compatible with all 2D tube formats
- Can be integrated into robotics
- > Detects missing tubes
- Cryoprotection function Active
 Data export in csv, html, Excel, Word,..
- Compatible with Windows 7, 8 and 10
- USB connection
- Dimensions W x D x H : 97 x 137 x 160 mm



eppendorf

greiner

Cat. No.	Description	€ Excl. VAT
830370	Perception Scanner for tubes on rack	NC -

Eppendorf 2D barcode scanner

2D barcode scanner

- >2D barcoded tube readings in SBS rack format, or cryobox
- Cryoprotection: special coating to prevent condensation from cryotubes from interrupting reading
- Small footprint
- PC connection via USB
- ▶ Compatible : Windows Vista, Windows 7, Windows 8
- Export of data in TXT, JSON, XML, image format
- > Documentation (certificate, etc.) on the tubes available via the Eppendorf website based on the cryotube codes
- > Automated data transfer to sample management software (e.g. Eppendorf eLABInventory)
- eLABJournal available via the Eppendorf website: intuitive and flexible solution for documenting research and tracking research data as well as workflow processes in which the ULT freezer is used

Cat.No	Description	€ Excl. VAT
934475	Rackscan 2D barcode scanner with Rackscan analyzer software	NC -

934476

1D barcode scanner

ID barcode scanner for lateral reading of SBS racks and cryoboxes
Can be mounted on the 2D barcode scanner Cat. number 934475

Cat. No.	Description	€ Excl. VAT
934476	1D barcode scanner	NC -

Bio-One Cryo.S[™] barcode scanner

- Decoding time: 8 to 10 seconds per rack
- >Anti-condensation treatment: for reading frozen tubes and racks
- > Ability to enable or disable internal temperature control during exposure to frozen samples
- Supplied with software for reading a tube or rack of tubes with 2D barcode, power cable, USB cable, instructions for use
- > One touch function for barcode reading activation

Cat. No.	Description	Dimensions L x W x H (mm)	€ Excl. VAT
849070	Cryo.S™ barcode scanner	310 x 195 x 70	NC -

