

Socorex micropipettes

Model	Acura® Dilute 810	Acura® manual 825	Acura® manual 835	Acura® manual 826XS	Acura® manual 855	Acura® Electro 926, 936 et 956
Operating mode	Manual	Manual	Manual	Manual	Manual	Electronic
Highlights	Straw tip, two volume dispensing for rapid sequential dilutions	Lightweight, robust, easy care	Lightweight, robust, easy care, optional pasteur pipette and straw adapter	Short, very light, robust, easy care	Lightweight, robust, 360° rotation of the bottom part, easy care	Intuitive programming, six working modes, fast charging, long autonomy
Volume range	1 ml and 0.1 ml	0.1 µl - 1 ml	0.2 µl - 10 ml	0.1 µl - 1 ml	0,5 µl - 350 µl	0.1 - 10,000 µl
Available options	1-channel	1-channel	1-channel	1-channel	8 channel, 12 channel	27 adaptable volumetric modules
Autoclavable Parts	Entirely	Entirely	Entirely	Entirely	Entirely	Lower part

Micropipette Acura® Dilute 810



- ▶ Easy sampling in narrow containers thanks to the straw tip
- ▶ Dispensing without adjustment of 1 and 0.1 ml from a single pipetting with two pre-calibrated volumes
- ▶ Allows easy sequential dilution to 1:10
- ▶ Independent calibration of each volume
- ▶ Ergonomic shape, effortless pipetting
- ▶ Adjustable ejector
- ▶ Interchangeable PE protective filter
- ▶ Easy care and clean, autoclavable without disassembly

Cat. No.	Description	Units/bag	Units/carton	€ Excl. VAT carton
041120	Fixed volume sampling pipette Acura® 810	-	1	NC -
041110	Sterile straw pipette 1,1 ml	25	1000	NC -
041115	Sterile straw pipette 2,2 ml	25	1000	NC -
061968	Replaceable PE protective filter	-	3600	NC -

Micropipettes Acura® manual 825, 826XS, 855 and 835



- ▶ Lightweight, optimised grip and smooth piston activation
- ▶ Volume control with precision stop
- ▶ Rotating cap prevents accidental volume changes
- ▶ Fast and safe calibration by the user
- ▶ Autoclavable without disassembly
- ▶ Robust, shock and UV resistant
- ▶ Warranty: 3 years

Micropipette Acura® 825 and 826XS

- ▶ Acura® 826XS : shorter and lighter with a thinner bottom section to facilitate pipetting into microtubes

Volume (µl)	Increment (µl)	Suitable tips	Acura® 825		Acura® 826XS	
			Cat.No	€ Excl. VAT	Cat.No	€ Excl. VAT
0,1 - 2	0.002	10 µl	061104	NC -	061700	NC -
0,5 - 10	0.01	10 µl	061105	NC -	061701	NC -
1 - 10	0.01	200 µl	061106	NC -	061702	NC -
2 - 20	0.02	200 µl	061107	NC -	061703	NC -
5 - 50	0.1	200 µl	061108	NC -	061704	NC -
10 - 100	0.1	200 µl	061109	NC -	061705	NC -
20 - 200	0.2	200 and 350 µl	061110	NC -	061706	NC -
100 - 1,000	1	1000 µl	061111	NC -	061707	NC -