


Multi One Touch bi-component non-filter and filter tips

- > Two-component tips: upper part in soft elastomer
- > Optimal tightness and precision: ideal tips for multi-channel pipetting
- > Reduced muscle fatigue for the user



One Touch bi-component tips

Tip	Plastic rack		
	Cat.No	Units/box	€ Excl. VAT/box
20 µl non sterile	010300	10 x 96	NC -
20 µl sterile	010310	10 x 96	NC -
330 µl non sterile	010330	10 x 96	NC -
330 µl sterile	010340	10 x 96	NC -
1000 µl non sterile	010360	8 x 96	NC -
1000 µl sterile	010370	8 x 96	NC -

One Touch bi-component filter tips

Volume (µl)	Plastic Cell-Rack		
	Cat.No	Units/carton	€ Excl. VAT/carton
1 - 20	010320	10 x 96	NC -
10 - 100	010510	10 x 96	NC -
20 - 350	010350	10 x 96	NC -
100 - 1250	010380	8 x 96	NC -

Multi Tips

> A guarantee of quality



> Advantages

- A tip especially designed for research
- A graSoft Tip: made of a softer polypropylene, these tips make it easier to attach to the micropipette and reduce the effort.
- A wide range of packaging and presentation options - see below.

> Standard tip or Low Binding tip

Multi tips offer two surface qualities:

- the "standard" surface in virgin polypropylene;
- the "Low Binding" surface: the tip, made hydrophobic, reduces surface tensions. This property has two advantages: maximum sample recovery and much lower DNA and protein retention than with an ordinary tip.

All Cat. No. : Next pages

Multi Tips, bulk



In a bag with a zip to guarantee re-sealing.

Hinged rack



Rack with hinged lid, hinge and lock.

CellRack Racked tips



- > Rack with hinged or covering lid
- > Individual closed cells for each tip
- > Design eliminates the risk of tip-to-tip contamination
- > Easy insertion of the tips on the micropipette
- > Enhanced rack stability
- > Autoclavable

Refill tips on Platinum tray

- > Sturdy pyramid hinged rack, stackable, repeatedly autoclavable
- > Saves space and waste
- > Transparent packaging to allow visualization of the tips



Refill Tips on QuickRack



- > 96 tips on a tray
- > Simply take a tray and place it on any empty tip rack to have a loaded rack
- > Fast, economical, small footprint
- > Allows reuse of empty Cell Rack

TipStation racked Tips



- > 10 ready-to-use tip racks in a minimum of space
- > Simply remove one level after using the tips
- > High stability bottom
- > Ideal for multi-channel micropipettes
- > Two-compartment lid