

Incubator-Clean™



PanReac
AppliChem
ITW Reagents

Disinfectant solution for incubators, sterile cell culture benches and molecular biology laboratories. Prevents contamination and growth of fungi and spores, yeast, bacteria and spores, mycoplasma and viruses (including HIV and Hepatitis B).

- **Active ingredient:** benzylammonium quaternary
- **Nontoxic and biodegradable**
- **Compatible with most work surfaces:** aluminium, chrome, nickel steel, polycarbonate..
- **Sterile, filtered solution in a spray bottle**
- **Storage:** at room temperature
- **Stability:** minimum 3 years
- **Recommendations for use :**
 - incubators: spray incubators every 15 days
 - dry without leaving any residue, so there is no need to rinse
 - sterile benches: to be used before each manipulation on the bench

Cat. No.	Type	€
670931	Incubator Clean 500 ml	NC -

Incuwater-Clean™



Disinfectant for use in incubator water tanks.

- **Non-toxic, non-volatile**
- **Does not damage stainless steel or copper trays in CO₂ incubators**
- **Conditions of use:** add 50 ml of Incuwater-Clean to the water tank with sterile water - repeat every fortnight

PanReac
AppliChem
ITW Reagents

Cat. No.	Type	€
670934	Incuwater Clean 100 ml	NC -

Aquabator-Clean™ (100x)



Disinfectant for use in water baths against bacteria and fungi.

- **Operating conditions :** 10 ml per litre of water - change water every 4-6 weeks and add aquabator-clean every 1-2 weeks
- **Stable at 100 °C**
- **Biodegradable**

PanReac
AppliChem
ITW Reagents

Cat. No.	Type	€
670935	Aquabator Clean (100x) 250 ml	NC -

Disinfectant for water baths Réa-Stabil



Concentrated bactericidal disinfectant for water baths and thermostatic baths. The results are even better with distilled or demineralised water.

- **Efficiency :** 2 to 3 months
- **Water colouring:** blue
- **Very economical:** 1 ml per litre (10 ml for a 10 litre bath)
- **Biodegradable**
- **Components:** biguanides surfactants and dye

Cat. No.	Description	€
074022	Bottle of 200 ml Réa-Stabil	NC -

Autoclave cleaner



- **Weak acid to be used diluted to clean autoclaves**
- **For use on walls and also for cleaning water tanks**

Cat. No.	Description	€
032204	Autoclave cleaner 1 l	NC -

Mycoplasma - prevention



Biogard S

- **Surface disinfectant**
- **Compatible with all types of surfaces**
- **Free of Alcohol, phenol and formaldehyde**
- **Non-toxic and biodegradable**

Cat. No.	Description	Vol.	€
500176	Biogard S Surface disinfectant spray	1 l	NC -



Biogard A disinfectant for CO₂ incubator water tank

- **Non-toxic, non-volatile**
- **Compatible with stainless steel**

Cat. No.	Description	Vol.	€
500177	Biogard A disinfectant for waterbaths and CO ₂ incubators	100 ml	NC -