

Corning® CellSTACK® culture chambers

- Five sizes: from 1-stack (636 cm²) to 40-stacks (25 440 cm²)
- Several surface treatments available: standard cell culture treatment, Corning CellBIND® surface for improved cell adhesion, or ULA (Ultra-Low Attachment) surface for minimal cell adhesion
- Sterile
- Non-pyrogenic
- Manufactured in the USA under GMP conditions
- Easy to use
- Wide opening with filter cap
- Format can be automated
- **Corning CellBIND® surface**
 - Ideal for reducing the addition of serum
 - Better attachment for better cell performance
- **ULA (Ultra-Low Attachment) surface: see chapter 3D**



Cat. No.	Surface area (cm ²)	Number of trays	Treatment	Units/bag	Units/carton	€ Excl. VAT/carton
003268	636	1	TC	1	8	NC -
003330	636	1	CellBIND®	1	8	NC -
003269	1272	2	TC	1	5	NC -
003310	1272	2	CellBIND®	1	5	NC -
003313	3180	5	TC	1	8	NC -
003319	3180	5	TC	1	2	NC -
003311	3180	5	CellBIND®	1	2	NC -
003270	6360	10	TC	1	2	NC -
003271	6360	10	TC	1	6	NC -
003312	6360	10	CellBIND®	1	2	NC -
003320	6360	10	CellBIND®	1	6	NC -
003272	25440	40	TC	1	2	NC -
003321	25440	40	CellBIND®	1	8	NC -

