

## Nunc™ LabteK II™ Culture chambers

### On slide, separable with adhesive




- Surface of the slide culture-treated
- 2 versions: culture-treated soda-lime glass or CC<sup>2</sup> glass, chemically modified to be positively charged reproducing the properties of a poly-D-lysine coating
- Dimensions L x W : 75 x 25 mm
- Superfrost™ marking area on the soda-lime glass slide and light blue marking area on the CC<sup>2</sup> glass slide
- Minimum fluorescence
- Comes with removal tool that separates the chamber from the slide
- The biocompatible adhesive acts as a barrier between the chamber and the slide when the two are assembled; once the chamber is removed, the slide has a hydrophobic edge that isolates the wells

Cat. No.	Number of compartments	Slide material	Units/ carton	€ Excl. VAT/ carton
055076	1	Soda-lime glass	16	NC -
055077	2	Soda-lime glass	16	NC -
055078	4	Soda-lime glass	16	NC -
055079	8	Soda-lime glass	16	NC -
055596	1	Glass CC <sup>2</sup>	16	NC -
055597	2	Glass CC <sup>2</sup>	16	NC -
055598	4	Glass CC <sup>2</sup>	16	NC -
055599	8	Glass CC <sup>2</sup>	16	NC -

- Polystyrene chamber glued to a microscope slide or coverslip
- Biocompatible adhesive, hydrophobic edges for each compartment
- Lid with handle on top
- Recommended volume depending on the surface

Number of compartments	Culture area (cm <sup>2</sup> )	Working volume (ml)/ chamber
1	8.6	2 - 4.5
2	4	4386 - 2
4	1.7	0.5 - 1
8	0.7	0.2 - 0.5

### On coverslip, non-separable with adhesive



- Borosilicate glass 1.5
- Dimensions L x W : 55 x 24 mm, thickness: n°1,5 (0,016 - 0,019 mm)
- The chamber cannot be separated from the coverslip
- Optimal for high power inverted microscopy observations
- Supplied on trays suitable for incubator use and stackable

Cat. No.	Number of compartments	Slide material	Units/ carton	€ Excl. VAT/ carton
055084	1	Borosilicate glass	16	NC -
055085	2	Borosilicate glass	16	NC -
055086	4	Borosilicate glass	16	NC -
055087	8	Borosilicate glass	16	NC -

## Merck Millipore® Millicell EZ slide



- Polystyrene chamber glued to a glass microscope slide
- The chamber can be easily separated from the slide, without tools
- Numbered wells
- Recommended volume depending on the surface

Number of compartments	Culture area (cm <sup>2</sup> )	Working volume (ml)/ chamber
4	1.7	0.5 - 1.7
8	0.7	0.2 - 0.7

Cat. No.	Number of compartments	Slide material	Units/ bag	Units/ carton	€ Excl. VAT/ carton
051988	4	Glass	4	16	NC -
051989	4	Glass	16	96	NC -
051990	8	Glass	4	16	NC -
051991	8	Glass	16	96	NC -