

CryoGen® Storage Box with Barcode

ClearLine® Polycarbonate boxes for cryotubes storage

- Polycarbonate box, grid and lid
- Transparent cover, with digital marking of each slot
- Available for all ClearLine® cryotube volumes, in :
  - to accommodate 25 and 81 tubes: coloured grid, white bottom
  - to accommodate 100 tubes: white bottom with pins
- Can be used from -196 °C to +121°C, autoclavables 20 min
- Do not store in the liquid phase of liquid nitrogen
- Laser engraved barcodes



Cryogen Box 2D bottom: 2D barcode and openings\*


	For	25 cryotubes 1-2 ml	81 cryotubes 1-2 ml	81 cryotubes 3-4-5 ml	100 cryotubes 2D - 2 ml	100 cryotubes 2D - 5 ml
	Dimensions (mm)	76 x 76 x 52	132 x 132 x 52	132 x 132 x 94	132 x 132 x 52	132 x 132 x 94
<b>Cryogen® Box 2D :</b> Standard barcode + 2D Datamatrix code	White bottom/ blue grid	-	390717	390721*	-	-
	White bottom/ yellow grid	-	390718	390722*	-	-
	White bottom/ red grid	-	390719	390723*	-	-
	White bottom/ green grid	-	390720	390724*	-	-
	White bottom with pins	-	-	-	390820**	390821**
	€ Excl. VAT	-	NC -	NC -	NC -	NC -
<b>Cryogen® Box 1D :</b> Standard Barcode	White bottom/ blue grid	390804	390808	390812*	-	-
	White bottom/ yellow grid	390805	390809	390813*	-	-
	White bottom/ red grid	390806	390810	390814*	-	-
	White bottom/ green grid	390807	390811	390815*	-	-
	White bottom with pins	-	-	-	390816**	-
	€ Excl. VAT	NC -	NC -	NC -	NC -	-

\* The CryoGen® Box for 81 tubes 3 to 5 ml is not compatible with 4 ml cryotubes with internal screw cap.

\*\* With internal screw cap only.

Numbering grids for Cryogen® Box 1D

- In polycarbonate
- To be placed in the bottom of Cryogen Box 1D 25 or 81 places
- Reflect the numbering printed on the lid of the boxes
- Minimise placement errors



Cat. No.	Description	€ Excl. VAT/ 10 units
390817	Numbered plate 1 to 25 for 25 places boxes	NC -
390818	Numbered plate 1 to 81 for 81 places boxes	NC -

STORAGE  
SCREENING


ClearLine® 2D barcode cryotubes

ClearLine® Quality

- » Made of medical grade virgin polypropylene
- » Can be used down to -196 °C
- » ATP-free
- » USP Class VI (non-cytotoxic and non-haemolytic)
- » Sterilised by Beta rays, SAL 10<sup>-6</sup>



- Graduated with inscription surface
- Autoclavable in steam at 121°C for 20 minutes.



- 2D Datamatrix barcode on the bottom of the tube, matching the printed barcode
- Pre-printed barcode on every tube (Code 128)
- Readable both ways, by all scanner brands



- Can be used with capper/decapper.
- New bi-injected screw cap (polyethylene and silicone)
- Can be screwed and unscrewed with one hand
- The bottom of the tubes can be twist-locked into the base of CryoGen® Box 2D and in the polycarbonate racks



	External screw		Internal screw	
Volume (ml)	2	5	2	5
H x ø (mm)	45.6 x 12.6	87.6 x 12.6	48.2 x 12.6	90.3 x 12.6
Cat. No.	390713	390714	390715	390716
€ Excl. VAT, 50 units	NC -	NC -	NC -	NC -