

PIPETMAN® M Connected





- Motorised micropipette requiring almost no pipetting force and preserving the gestures associated with a mechanical pipette
- ➤ Lightweight, low ejection force and ergonomic design
- **▶**Lockable volume
- >Simple and intuitive operation with digital display on OLED screen
- >4 pipetting modes: pipette, repeat, mix and invert
- ➤ The PIPETMAN[®] M Connected software allows you to create and transfer up to 10 custom protocols
- ➤ Rechargeable battery in just 3 hours for up to 900 pipetting cycles
- ➤ Can be used while charging
- ➤ Can be connected to the optional Trackman® Connected to help reduce human error, track pipetting or archive calibration data
- ➤ Autoclavable lower part and ejector
- ➤ Volumetric specification according to ISO 8655-2
- ➤Warranty: 1 year

Single channel


| Volume (µl) | Suitable tips | Cat. No. | € Excl. VAT |
|----------------|---------------------|----------|-------------|
| 0.5 - 10 | D(F)(L)10 | 067516 | NC - |
| 2 - 20 | D200, DF30 | 067517 | NC - |
| 5-100 | D200,DF100 | 067518 | NC - |
| 20 - 200 | D(E) D(E) | 067519 | NC - |
| 20 - 300 | D(F)200, D(F)300 | 067520 | NC - |
| 100 - 1,200 | D(F)1,000, D(F)1200 | 067521 | NC - |
| 500- 5,000 | D5,000 | 067522 | NC - |
| 100- 10,000 | D10 ml | 067523 | NC - |

Multichannel

| | | 8-channel | | 12-channel | |
|----------------|---------------------|-----------|----------------|------------|----------------|
| Volume (µl) | Suitable tips | Cat. No. | € Excl. VAT | Cat. No. | € Excl. VAT |
| 0.5 - 10 | D(F)(L)10 | 067524 | NC - | 067525 | NC - |
| 1 - 20 | D(F)L10, D200, DF30 | 067526 | NC - | 067527 | NC - |
| 10-100 | D200,DF100 | 067528 | NC - | 067529 | NC - |
| 20 - 200 | D(F)200, D(F)300 | 067530 | NC - | 067531 | NC - |
| 10 - 300 | | 067532 | NC - | 067533 | NC - |
| 50 - 1,200 | D(F)1,000, D(F)1200 | 067534 | NC - | 067535 | NC - |

Trackman® Connected





- ➤ Tablet with exclusive accessories and applications to connect to the pipetman® M connected
- Tracks pipetting on 96 and 384 well plates to avoid errors
- >Interacts in real time with the PIPETMAN® M connected to guide the user through the protocol
- ➤ Recording and sharing of results and protocols
- ➤ Monitoring of usage statistics and calibration data history
- Supplied with a microplate holder, microplate adapter and PCR strip; and a sensor to monitor temperature, pressure and humidity during handling

| Cat. No. | Description | € Excl. VAT |
|----------|---------------------|----------------|
| 400115 | Trackman® Connected | NC - |