

ART filter tips



All brands of micropipettes have a perfectly fitting ART tip. The Softfit tip has a reinforced seal. The hydrophobic filter prevents any accidental backflow of liquid into the filter and pipette.

All ART tips are sterile, (except those in bulk packaging). The ART tip is available in a standard version, or a low retention version for improved sample recovery.

A wide range of presentations:

- on a single rack with a cover lid;
- on individual hinged-lid racks;
- In bulk ;
- in individual refills: the sterile tips are placed on a tray under an individual blister pack;
- in tray refills, ready to use

- 1) Hold the top of the blister pack and remove the peel-off sheet.
- 2) Place the blister pack on the rack.
- 3) Press.

Sterile refills



Model	Vol. (µl)	Length (mm)	Features	Compatibility	Pack unit
ART 10 Reach	0 – 10	43.3	Ideal for micro-volumes, long length suitable for microtubes	Gilson P2-P10, Finnpiquette 2 and 10 µl micro, Sartorius Biohitfamily 10 µl, Eppendorf 2.5 µl	10 x 96
ART 10	0 – 10	31.3	Low volume	Gilson P2-P10, Finnpiquette 2 and 10 µl, Sartorius Biohitfamily 10 µl	10 x 96
ART 20E	0 – 20	47.7	Low volume	Eppendorf 10 µl, Socorex 10 µl	10 x 96
ART 10F	0 – 10	58.8	Low volume	Finnpiquette 10 µl universal, Sartorius Biohitfamily 10 µl	10 x 96
ART 20P	0 – 20	50.7	Universal	Gilson P20, Eppendorf 20 µl, Socorex 20 µl	10 x 96
ART 20	0 – 20	45.7	-	Eppendorf 10 µl	10 x 96
ART 50U	0 – 50	52	Softfit, with reinforced rack, special multi-channel	Multichannel (Labsystems 5 - 50 µl, Sartorius Biohitfamily 5 - 50 µl, Eppendorf 5 - 50 µl)	10 x 96
ART 100	1 – 100	58.8	-	Finnpiquette 40 µl Multichannel (Labsystems 5 - 50 µl, Sartorius Biohitfamily 5 - 50 µl)	10 x 96
ART 100E	1 – 100	50.7	Universal	Gilson P100, Eppendorf 100 µl, Sartorius Biohitfamily 50 µl and electronic 100 µl	10 x 96
ART 200	1 – 200	58.8	Universal	Gilson P200, Eppendorf 200 µl, Finnpiquette 200 µl, Sartorius Biohitfamily 200 µl	10 x 96
ART 200 REACH	1 – 200	90.9	Long length	Gilson P200, Eppendorf 200 µl, Finnpiquette 200 µl, Sartorius Biohitfamily 200 µl	8 x 96
ART 200U	1 – 200	52	Softfit, with reinforced rack, special multi-channel	Multichannel (Labsystems 50 - 300 µl, Eppendorf 50 - 300 µl), any 200 µl pipette	10 x 96
ART 200G	5 – 200	49.7	Large orifice for fragile cells, DNA and viscous solutions	Gilson P200, Eppendorf 200 µl, Finnpiquette 200 µl, Sartorius Biohitfamily 200 µl	10 x 96
ART 300	1 – 300	58.4	Multi-channel special	Multichannel (Sartorius Biohitfamily 250 µl, Finnpiquette 300 µl), Eppendorf Research Pro 300 µl	10 x 96
ART 500	50 – 500	79	Softfit	Eppendorf 50 - 250 µl, Sartorius Biohitfamily 50 - 500 µl	8 x 100
ART 1000	100 – 1000	79	Softfit	Gilson P1000, Eppendorf 1000 µl, Finnpiquette 1000 µl, Sartorius Biohitfamily 1000 µl	8 x 100
ART 1000E	100 – 1000	102.2	Long length	Gilson P1000, Finnpiquette 1000 µl, Sartorius Biohitfamily 1000 µl, Eppendorf 1000 µl	8 x 100
ART 1000G	100 – 1000	85.4	Large orifice for fragile cells, DNA and viscous solutions	Gilson P1000, Eppendorf 1000 µl, Finnpiquette 1000 µl, Sartorius Biohitfamily 1000 µl	8 x 100
ART 1000 REACH	100 – 1000	118	Long length, Softfit	Gilson P1000, Eppendorf 1000 µl, Finnpiquette 1000 µl, Sartorius Biohitfamily 1000 µl	8 x 96
ART 1250 REACH	100 – 1250	102.2	Long length	-	8 x 96
ART 5000	1000 – 5000	117	-	Gilson P5000	Sterile bag of 250 tips
ART GEL 20P	0.5 – 20	70.9	A flexible capillary for gel coating	Gilson P20, Eppendorf 20 µl, Finnpiquette 10 µl	10 x 96
ART GEL	0.5 – 100	70.9	A flexible capillary for gel coating	Gilson P100, Eppendorf 100 µl, Finnpiquette 40 µl	10 x 96
ART XLP	5 – 100	91.2	Long length; reaches the bottom of 12 x 75 mm tubes	Gilson P100, Eppendorf 100 µl, Finnpiquette 40 µl	8 x 96
ART XLP 200	5 – 200	91.2	Long length; reaches the bottom of 12 x 75 mm tubes	Gilson P200, Eppendorf 200 µl, Socorex 1000 µl	8 x 96
ART XLG	5 – 200	91.2	Long length and wide opening for cells, DNA and viscous solutions	Gilson P200, Eppendorf 100 µl, Eppendorf 200 µl, Finnpiquette 40 µl	8 x 96
ART 1000 XL	100 - 1000	89	Universal	Gilson P1000, Eppendorf 1000 µl, Sartorius Biohitfamily 1000 µl Finnpiquette 1000 µl	8 x 96

“Essential” Kits

- “Essential” hinged rack kits contains:
- ART 200 (Cat. No. 038019) x 4
 - ART 10 reach (Cat. No. 038001) x 4
 - ART 1000XL (Cat. No. 041539) x 4

- “Essential” refill insert kits contains :
- ART 200 (Cat. No. 038019) x 4
 - ART 10 reach (Cat. No. 038001) x 4
 - ART 1000XL (Cat. No. 041539) x 4

Cat.No	Description	€ Excl. VAT
041604	"essential" ART tips Kit: 4 hinged racks of each tip Cat. No. 038019, Cat. No. 038001, Cat. No. 041539	NC -
041605	"essential" ART tips Kit: 4 refill inserts of each tip Cat. No. 038019, Cat. No. 038001, Cat. No. 041539	NC -