

## Movebox - transport solution for biological samples

### Movebox

- Complies with ADR 2009 - instruction P650 for the transport of category B biological samples
- Complies with GBEA II 2002 and NF EN ISO 15189 2007
- The system can hold up to 77 tubes of 5, 7 or 10 ml
- Rack with 2 integrated gripping handles
- Ergonomic solution

### Movebag cooler bag

- Capacity of up to 4 Moveboxes, i.e. 308 tubes
- Internal resistant pocket for eutectic plates
- Easy to decontaminate



### Eutectic plates for Movebag

- Optimised temperature management: ambient (+22°C) or +4°C

### Tube rack

- For tubes of 13 or 16 mm diameter
- 2 handles for easy handling of the rack
- Easy decontamination

### Absorbent for Movebox box

- Sheet with gel micro-particles to absorb liquid in the event of a leak

Cat. No.	Description	€ Excl. VAT
368700	Movebox kit (Movebag cooler bag, 4 Movebox boxes with holder for 13 mm diameter tubes and 2 absorbents/box)	NC -
368701	Movebag cooler bag for 4 Movebox boxes	NC -
368702	Movebox box with holder for 13 mm diameter tubes and 2 absorbents	NC -
368703	Movebox box with holder for 16 mm diameter tubes and 2 absorbents	NC -
368704	Box Movebox, Natural	NC -
366907	Box Movebox, White	NC -
361575	Box Movebox, Yellow	NC -
361577	Box Movebox, Purple	NC -
361578	Box Movebox, Green	NC -
361579	Box Movebox, Pink	NC -
368705	Rack for 13 mm diameter tubes	NC -
368706	Rack for 16 mm diameter tubes	NC -
368707	Eutectic plate 22 °C, 200 g	NC -
368708	Eutectic plate 22 °C, 1100 g	NC -
368709	Eutectic plate 4 °C, 200 g	NC -
368710	Eutectic plate 4 °C, 1100 g	NC -
368711	Absorbent for Movebox box	NC -

## Freez'Pack ice block



- To keep products, liquids or medicines cool
- Safe (non-toxic refrigerant)
- Reusable at will
- Leave in the freezer for 2-4 hours before use
- Produces cold as long as it is frozen
- Place the Freez'Pack on top of the products to be kept refrigerated

Cat. No.	Model	For Coolers, volume	Dimensions (mm)	Units/box	€ Excl. VAT/box
071110	M5	1 to 6 l	150 x 76 x 24	2	NC -
071111	M10	6 to 10 l	180 x 95 x 32	1	NC -
071112	M 20	10 to 19 l	200 x 170 x 30	1	NC -
071113	M 30	19 to 32 l	273 x 200 x 30	1	NC -

## New Ice cushion



Ice cube in the form of a cushion filled with ultra-pure water. Ideal for storing chilled samples without the inconvenience of melting ice cubes.

- Formed from a high purity FDA approved polyethylene bag
- Filled with ultra pure water
- Reusable, washable and refreezable

Cat. No.	Description	Units/box	€ Excl. VAT/box
053096	New Ice, Ice cushion	45	NC -