

## Nalgene™ Singles use Top-Filters



### Features PES FASTER

- PES FASTER membranes : Supor Mach V membrane
- Rapid-flow membrane: new manufacturing technology for faster and more efficient filtration
- Asymmetrical
- Low protein binding rate



Fits onto ISO, Pyrex or Schott glass bottles with a 33 or 45 mm neck. PES (0.45 µm) or PES FASTER (0.2 µm) membranes, and SFCA membranes.

Membrane	Volume (ml)	Porosity (µm)	ø (mm)	33 mm neck		45 mm neck	
				Cat. No.	€ 12 units	Cat. No.	€ 12 units
PES FASTER	150	0.2	50	029670	NC -	029671	NC -
PES	150	0.45	50	028599	NC -	028609	NC -
PES FASTER	500	0.2	75	029668	NC -	029669	NC -
PES	500	0.45	75	028601	NC -	028611	NC -
PES FASTER	1000	0.2	90	029672	NC -	029673	NC -
SFCA	150	0.2	50	028603	NC -	028613	NC -
SFCA	150	0.45	50	028604	NC -	028614	NC -
SFCA	500	0.2	75	028605	NC -	028615	NC -
SFCA	500	0.45	75	028606	NC -	028616	NC -
SFCA	1000	0.2	90	028607	NC -	028617	NC -

## Receiving bottles for Nalgene™ filtration units



Polystyrene bottles with polyethylene screw cap for Nalgene™ filtration units.

### ➤ Sterile ➤ 45 mm neck

Cat. No.	Volume (ml)	Units/box	€/box
029190	150	24	NC -
029191	250	24	NC -
029192	500	12	NC -
029193	1000	12	NC -

## Nalgene™ Analytical Filtration units



Polypropylene filter unit with removable membrane and receiving flask or funnel.

- Cellulose nitrate membrane 0.45 µm or 0.2 µm
- The membrane is easily removed
- Volume : 150 ml
- Gamma sterilised, single use

Volume (ml)	150	150	150
Receiver	Flask	Flask	funnel
Membrane	NC	NC	NC
Porosity	0.2 µm	0.45 µm	0.45 µm
Colour	White	White with grid	White with grid
Units/bag	12	12	12
Units/box	72	72	72
Cat. No.	028032	029335	029336
€ the box	NC -	NC -	NC -

## Nalgene™ Analytical Filtration Funnels



Economical, graduated, with a membrane for filtration and microbiological analysis. The top polypropylene part of the funnel can be removed for immediate access to the membrane.

- Cellulose nitrate membrane, approved for water control
- Sterile, single use
- White or grey membrane (turning black when wet)

Cat. No.	Volume (ml)	Porosity (µm)	Membrane colours	Gridded	Units/box	€/box
029378	100	0,2	White	Without	50	NC -
028414	100	0.45	White	Grey	50	NC -
028004	100	0.45	Grey	White	50	NC -
029379	250	0,2	White	Without	50	NC -
028415	250	0.45	White	Gray	50	NC -

### Accessories

Cat. No.	Description	€
029380	Autoclavable polypropylene funnel rod only (25 units)	NC -