

## Multicaster




RunOne gel casting accessories for high density screening (96 well agarose gels). The Multicaster system consists of a tray with notches and combs designed for casting agarose gels in a few minutes without seals or adhesives. The trays are UV transparent. 4 models available.

### Multicaster

Type	8 rows	12 rows	Preparatory gels	HLA Multicasting
L x W (mm)	116 x 118	116 x 118	116 x 118	116 x 89
Possible deposit lines	8 x (12 + 1) / 4 x (24 + 1)	12 x (8 + 1) / 6 x (16 + 2)	2	8
Recommended pipette	Single channel / 12-channel	8-channel	Single channel	Single channel / 12-channel
Deposit volume (µl)	30 – 35 / 10 – 15	30 – 35 / 10 – 15	2 x 250 + marker / 3 x 75 and 2 x 115	30 – 35 / 10 – 15
Number of notches for combs	8	12	2	8
Number of combs included (1 mm thick)	4	6	2	8
Reversible double comb types included	12 + 1 / 24 + 1	8 + 1 / 16 + 1	2 / 5 + markers	12+1 / 24+1
Maximum number of deposits	192 / markers	192 / markers	2 / 5 + markers	192 / markers
Colour code	Turquoise	Blue	Orange	Yellow
Cat. No. complete unit	038415	038416	038417	039080* / 038418**
€ Excl. VAT.	NC -	NC -	NC -	NC - / NC -
Cat. No. Multicaster + RunOne Kit	035459	035460	035461	035462B***
€ Excl. VAT.	NC -	NC -	NC -	NC -
Cat. No. replacement tray	035468	035468	035469	035469
€ Excl. VAT. 2 units	NC -	NC -	NC -	NC -
Cat. No. comb	038449	039033	039070	039071
€ Excl. VAT. 2 units	NC -	NC -	NC -	NC -

\* Supplied with 2 combs \*\* Supplied with 8 combs \*\*\* Supplied with the RunOne 6-compartment casting system and a red viewing plate

## Reversible combs for RunOne



### T-Blue comb:

To be used when preparing three gels, one in the large bowl, and the other two in the small ones.

### Reversible :

- one side 1.5 mm thick: 25/12/12 teeth ;
- one side 1 mm thick: 17/8/8 teeth

### T-Ambre Comb:

#### Preparation. Reversible :

- one side 1 mm thick 25/12/12 wells ;
- one side 1 mm thick, 10+2 wells, 2+2 preparatory wells, 1+2 preparatory wells.

Cat. No.	Description	€ Excl. VAT
039075	T-Blue Reversible comb (2 units)	NC -
039076	T-Ambre Reversible comb (2 units)	NC -
039077	T-Bleu + T-Ambre Reversible comb	NC -



Find RunOne compatible pre-cast gels on [www.dulis.nl](http://www.dulis.nl)