

Sealing Films and rolls

Cat.No	Photo	Model	Features	Packaging	For	€ Excl. VAT
016462F	A	Clear Seal	Optical Film for QPCR	500 m x 78 mm Roll	ALPS 300, 3000 and A4S	NC -
044770	A		Effective from -80 °C to 80 °C	350 m x 115 mm Roll	Velocity 11 PlateLoc	NC -
016427F	B		or 110 °C with pressure lid	100 foils	ALPS 50V or Sealer 4S3	NC -
016276F	A	Clear Weld Seal Mark II	Optical film for QPCR and storage	610 m x 78 mm Roll	ALPS 300, 3000 and A4S	NC -
044771	A		Effective from -80 °C to 110 °C	500 m x 115 mm Roll	Velocity 11 PlateLoc	NC -
044720	B		Good resistance to solvents (DMSO)	100 foils	ALPS 50V or Sealer 4S3	NC -
016274F	A	Clear Seal 3730	To be used with PE, PS and PP plates	610 m x 78 mm Roll	ALPS 300, 3000 and A4S	NC -
044773	A		Effective from -80 °C to 90 °C or 110 °C with pressure lid - pierceable	500 m x 115 mm Roll	Velocity 11 PlateLoc	NC -
044772	A		Recommended for the ABI 3730 sequencer	Roll of 1000 pre-cut foils	ALPS 50V or Sealer 4S3	NC -
016273F	C	Peel Seal	Peelable film for PP plates	600 m x 78 mm Roll	ALPS 300, 3000 and A4S	NC -
044767	C		Effective from -80 °C to 110 °C	500 m x 115 mm Roll	Velocity 11 PlateLoc	NC -
016279F	D		For storage and PCR	100 foils	ALPS 50V or Sealer 4S3	NC -
044777	C	DMSO Resistant Peel Seal	Peelable film for PP plates	500 m x 78 mm Roll	ALPS 300, 3000 and A4S	NC -
044778	C		Effective from -80°C to 40°C - for storage.	500 m x 115 mm Roll	Velocity 11 PlateLoc	NC -
044779	D		Very good resistance to solvents (DMSO)	Roll of 100 pre-cut foils	ALPS 50V or Sealer 4S3	NC -
016679F	E	Pierce Seal	Pierceable film for PP plates	610 m x 78 mm Roll	ALPS 300, 3000 and A4S	NC -
044768	E		Effective from -20 °C to 120 °C	500 m x 115 mm Roll	Velocity 11 PlateLoc	NC -
016678F	D		Very good resistance to solvents, resealing on itself possible for PCR and storage	100 foils	ALPS 50V or Sealer 4S3	NC -
016272F	E	Foil Seal	For PP or PS plates - Peelable and pierceable	610 m x 78 mm Roll	ALPS 300, 3000 and A4S	NC -
044769	E		Effective from -20 °C to 110 °C	500 m x 115 mm Roll	Velocity 11 PlateLoc	NC -
016278F	D		Good solvent resistance - Resealable on itself possible for PCR and storage	100 foils	ALPS 50V or Sealer 4S3	NC -
044774	F	Thermal Bond	Peelable, very high sealing efficiency	610 m x 78 mm Roll	ALPS 300, 3000 and A4S	NC -
016277F	D		High degree of sample protection -20 °C to 110 °C - PCR and long-term storage	100 foils	ALPS 50V or Sealer 4S3	NC -

Manual Heat Sealer ALPS 30 ABgene®

- Compact (dimensions : 220 x 250 x 230 mm)
- Lightweight (approx. 4 kg)
- Fast heating : LED indicating sealing temperature
- Ergonomic
- Includes: Deepwell plate adapter, 96-well PCR plate adapter and standard microplate adapter

Cat. No.	Description	€ Excl. VAT
017230	Heat sealer ALPS30	NC -



STORAGE
SCREENING

Semi-automatic heat sealer ALPS 50V

- Heat sealer for microplate film
- For low to medium flow applications
- Compact size
- Audible alarm at the end of the sealing process
- Temperature and time setting via touch pad

Welding T° range	125 °C to 200 °C
Operating T°	18 °C to 30 °C
Humidity range	20 to 80 %
Dimensions L x D x H (mm)	220 x 321 x 425
Weight (kg)	7.2
Cat. No.	446879
€ Excl. VAT	NC -

