

## Eyewash and saline or pH Neutral wall stations

### > Sterile 0.9% sodium chloride saline solution

- > For rinsing and removing foreign bodies, dust, splinters, solvents, oils
- > Rinse until all foreign bodies are removed

Estimated rinsing time (min.)		
Flask volume	Single eyelet	DUO Eyelet
200 ml	2	-
500 ml	5	2
1000 ml	10	5

### Bottles only

Eyelet	Flask volume (ml)	Cat. No. saline solution	€	Cat. No. pH Neutral	€
Simple	200	060030B	NC -	941012B	NC -
Simple	500	030025B	NC -	-	-
Simple	1000	060725	NC -	060728	NC -
DUO	500	060724	NC -	060727	NC -
DUO	1000	060726	NC -	-	-

### Wall-mounted stations and boxes

- > With saline solution and/or pH Neutralizing Solution
- > The polystyrene boxes are dustproof and impact resistant and can be hung on walls

Cat. No.	Description	Saline solution	pH Neutral	€
060729	Wall station (wall bracket & pictogram)	1 x 500 ml	-	NC -
060028	Wall station (stand & pictogram)	2 x 500 ml	-	NC -
060730	PH Neutral Wall Station	-	1 x 1000 ml	NC -
060731	Maxi Wall Station	2 x 1000 ml	-	NC -
060740	DUO wall station	1 x 1000 ml DUO	-	NC -
941008	Plum Combi Station	1 x 500 ml	1 x 200 ml	NC -
060741	Plum Combi Station DUO	1 x 1000 ml DUO	1 x 500 ml DUO	NC -
060732	Wall-mounted box	2 x 500 ml	-	NC -
060733	DUO box set	2 x 1000 ml DUO	-	NC -
060735	Plum Combi-Box DUO	1 x 500 ml DUO	1 x 500 ml DUO	NC -
060736	Plum Combi-Box DUO MAXI	1 x 1000 ml DUO	1 x 500 ml DUO	NC -
060744	PS empty box for 200 ml to 500 ml bottles	-	-	NC -

### > pH Neutral : 4.9 % sterile phosphate buffer for neutralisation :

- > For rinse and neutralisation of acid and alkali chemical splashes
- > Rinse with pH Neutral solution for about 2 minutes, then rinse with saline solution

### > One eye flasks :

- > For one eye rinsing at a time

### > DUO flask

- > For simultaneous rinsing of both eyes



### Quick Rinse eyewash ampoules

- > With saline solution and/or pH Neutralizing Solution
- > The polystyrene boxes are dustproof and impact resistant and can be hung on walls

Cat. No.	Description	Content	€
060737	Quick Rinse wall dispenser	10 ampoules	NC -
060738	Quick Rinse refill for dispensers	4 cases x 5 ampoules	NC -
060739	Quick Rinse Mini cases	5 cases x 4 ampoules	NC -

### Accessories

Cat. No.	Description	€
060742	Belt bag for 200 ml bottle (delivered empty)	NC -
060029	Belt bag with 200 ml saline solution bottle	NC -

## Portable eyewash bottle

- > 235 ml bottle (isotonic solution)
- > Contains a buffered solution: neutralising effect on acidic and basic splashes
- > Easy to carry, pocket size
- > Phosphate-free
- > Dust protected eyelet
- > Generous rinsing
- > Shelf life : 4.5 years

Cat. No.	Description	€
039056B	Portable eyewash bottle 235 ml	NC -

